

Understanding the factors affecting sensing efficiency in SPR biosensing

Introduction

Acute lymphoblastic leukemia (ALL) is a disease where immature white blood cells originated from the bone marrow become cancerous, inducing proliferation in the blood flow and subsequently to organs in addition to prevent the proper function of other blood cells. One of the main chemotherapeutic agents proven for its efficacy for ALL is the E. coli L-asparaginase (EcAII) as a biological therapeutic agent. ¹ The challenge with undergoing such treatment is that the patient can potentially develop silent inactivation of the biological chemotherapeutic agent by generating antibodies to neutralize EcAII and therefore reducing treatment efficiency.²

In order to detect the presence of anti-EcAII in patients, EcAII is employed in an immunoassay format. This method was previously used with the P4SPR to detect the anti-therapeutic antibody, anti-EcAII antibody, in undiluted serum of child patients.³ Biosensing with SPR has the advantage of minimizing sample preparation and reducing analysis time. However, the main challenge of developing the assay was to immobilize native EcAII antigen with reproducible sensor efficiency.

SPR is a powerful bioanalytical technique in assay optimization as it provides real-time monitoring of label-free targets. In this case, the immobilization density and the level of response to the anti- EcAll antibody were carefully monitored. This allowed researchers to understand the underpinning mechanism of the immobilization process and determine the optimal protocol to build a more efficient sensor for a clinically relevant immunoassay.

Immobilization challenges of native EcAII



Figure 1 – Comparison of immobilization approaches and resulting sensing efficiency and sensitivity.

In SPR assay development, the immobilization approach plays a crucial role in ensuring optimal and reproducible sensing efficiency. The immobilization chemistry may affect the efficacy of antibody recognition if the immobilized receptors have obstructed binding sites (Figure 1). This poses as a challenge for SPR sensing, or any other immunoassay, where the sensitivity depends on the optimal proteinprotein binding interaction. Another key challenge affecting the assay performance during the immobilization is maintaining the quaternary structure. Alterations to the quaternary structure upon immobilization could affect the SPR signal and produce incorrect information. Both the protein density on the surface and quaternary structure are in play for building an EcAll immunosensor.

¹ Swain, A. L.; Jaskolski, M.; Housset, D.; Rao, J. K. M.; Wlodawer, A. Crystal structure of Escherichia coli L-asparaginase, an enzyme used in cancer therapy. Proc. Natl. Acad. Sci. U.S.A. 1993, 90, 1474–1478.

² Shinnick, S. E.; Browning, M. L.; Koontz, S. E. Managing

hypersensitivity to asparaginase in pediatrics, adolescents, and young adults. J. Pediatr. Oncol. Nurs. 2013, 30, 63–77.

³ Aubé, A.; Charbonneau, D. M.; Pelletier, J. N.; Masson, J.-F. Response Monitoring of Acute Lymphoblastic Leukemia Patients Undergoing I-Asparaginase Therapy: Successes and Challenges Associated with Clinical Sample Analysis in Plasmonic Sensing. ACS Sens. 2016, 1, 1358–1365.

Application Note

Assay optimization of protein-protein interaction





The EcAll forms a tetramer with the dimerization of intimate homodimers to form a functional homotretrame. ⁴ The homotetrameric structure exhibiting four identical active sites formed by complementation of intimate homodimers. Although several antigenic determinants of EcAll have been identified, its antigenicity in an immobilized form remains unknown.

To better understand the contribution of these factors on patient immunogenicity to EcAll, a heterogeneous surface immobilization of native EcAll was compared to a homogeneous one. The surface immobilization of EcAll was performed by coordination of N- or C-



terminal histidine (His)-tags of EcAll variants, tetrameric as well as monomeric (Figure 1).

Immobilization of anti-EcAll

First, the SPR gold sensor was coated with AffiCoat, a peptide self-assembled monolayer (SAM), to reduce nonspecific binding and introduce –COOH groups. The terminal -COOH are available for covalent binding with lysine residues on native EcAII via amine coupling with EDC/NHS. For metal chelation of terminal His-tag, Co-NTA was further attached onto the AffiCoat. With both surface preparation, the EcAII was in contact with the surface for 20 minutes then washed off with buffer.

The native E. coli L-asparaginase II was commercially available via Kidrolase, EUSA Pharma. The design and production of the N- and C-terminally His-tagged versions of recombinant EcAII developed for immobilization and SPR sensing were done in-house. Three His-tagged versions of EcAII were tested: N21and N26-EcAII which harbor 20 and 25 additional amino acids and C8-EcAII which harbor 8 additional amino acids (Figure 2).

Immunosensor SPR performances

SPR biosensing experiments were performed using the benchtop P4SPR device. The prepared immunosensor, both native and recombined EcAII, was first passivated to reduce nonspecific binding later in the test. Blank human serum was injected and in contact with the sensor for 10 minutes. Then, human serum spiked with 2 concentrations of polyclonal rabbit anti-asparaginase antibodies was injected, and binding was monitored for 20 minutes.

⁴ Epp, O.; Steigemann, W.; Formanek, H.; Huber, R. Crystallographic evidence for the tetrameric subunit structure of

L-asparaginase from Escherichia coli. Eur. J. Biochem. 1971, 20, 432–437.

Application Note

Assay optimization of protein-protein interaction



	receptor (Ag) immobilization				antibody (Ab) detection ^{b}	
protein	$\Delta\lambda_{\text{SPR}}$ (nm)	surface density (pmol cm ⁻²)	activity ^{a} (U mg ⁻¹)	distance c.t.c. (nm)	$\Delta\lambda_{SPR}$ (nm)	sensing efficiency (Ab/Ag binding ratio)
Kidrolase	1 ± 0.7	0.3 ± 0.2	310 ± 260	30.5	2.5 ± 0.9	2.3 ± 0.8
N21-EcAII	13.5 ± 0.6	3.3 ± 0.2	87 ± 7	8.5	12.2 ± 1.3	0.9 ± 0.1
N26-EcAII	12 ± 2	2.9 ± 0.5	95 ± 10	9.1	11.3 ± 0.8	0.9 ± 0.1
C8-EcAII	9.7 ± 0.8	2.5 ± 0.2	160 ± 87	9.9	8.3 ± 0.7	0.9 ± 0.03
C8-EcAII m	4.6 ± 0.8	$4.7 \pm 0.8^{\circ} (1.2 \pm 0.8)$	0	14.4	3.2 ± 0.5	$0.2 \pm 0.03^{c} (0.7 \pm 0.12)$
EcAII _{crystal} d		1.9		11.4		
EcAIIm _{crystal} ^e		2.1		10.8		

^aSpecific activity of immobilized EcAII receptors (Ag) was monitored on 9×9 mm gold-coated glass slides. All other immobilization and detection measurements were performed in the P4SPR instrument using 20×12 mm gold-coated prisms. ^bDetection using $150 \,\mu g \, mL^{-1}$ anti-EcAII antibody (Ab). ^cDensity or binding ratio of the monomer; values comparable to 4 monomers (1 equiv tetramer) are in parentheses. ^dUnit cell dimensions for ECAII in the tetrameric form (PDB 3ECA). ^eUnit cell dimensions for ECAII in the monomeric form (PDB 1NNS).

Table 1 - Summary of immobilization and immunosensing performance of EcAII receptors.



Figure 3 – Immobilization of EcAII variants. SPR sensograms of surface immobilization of EcAII variants was followed upon the injection of 40 μ g of protein (0.1 mg mL⁻¹) for Kidrolase (black) or His-tagged EcAII: N21-EcAII (green), N26-EcAII (blue), and C8-EcAII (red). The monomeric variant is identified.

Table 1 summarizes the immobilization characteristics of the various EcAll variants on the SPR chip and their bioanalytical performance. The surface density calculation is shown in appendix.

Native EcAll and other variants were immobilized onto the SPR sensor chip for assessment of the impact of the mode of immobilization on the receptor surface coverage (Figure 3). The native EcAll was immobilized on the sensing surface in a random manner by covalent cross-linking of lysine residues with EDC/NHS chemistry whereas His-tagged EcAll variants were immobilized with metal coordination on the sensor surface. The surface coverage determined for oriented His-tagged EcAll was 5 to more than 10 fold greater than covalent cross-linked EcAll (Table 1). In addition, the surface coverage with His-tagged EcAII was highly reproducible with relative standard deviation (RSD) of 3 to 18% compared to 50 to 70% for native EcAll. The low surface coverage of the C8-EcAll (monomeric or tetrameric) in comparison to the N-terminal variants could be explained by the fact that the C-terminal tag (8 residues) is less accessible for coordination. Moreover, results indicated that the monomeric C8-EcAll did not re-assemble into a tetrameric form whereas the tetrameric EcAll variants did not change in structure upon immobilization as demonstrated by corresponding immunosensor SPR response, low SPR shift for the monomeric form and high SPR shift for tetrameric forms (Figure 1.).

Additionally, the P4SPR data showed that the Nterminal modified His-tagged EcAll had greater surface density relative to both Kidrolase and C8-EcAll which resulted in better detection of anti-EcAII detection as shown at 15 and 150 µg/mL. N26-EcAll and C8-EcAll demonstrated 2-fold increase in SPR detection than randomly oriented native EcAll. The signals were further boosted with oriented tetrameric N- and Cterminally tagged EcAll providing 5-fold and 3-fold greater detection signal than Kidrolase and oriented monomeric C8-EcAll respectively (Table 1). Regarding the sensing reproducibility, RSD of 35% of Kidrolase was also higher than the His-tagged variants with RSD of 5 to 20% at both target antibody concentrations (Table 1). Correlation was made between the immunosensing signal and the surface coverage of EcAll receptors with R² of 0.9571 at the higher antibody

Application Note Assay optimization of protein-protein interaction



concentration with the oriented N-terminally Histagged variant (Figure 4).



Figure 4 – Efficiency of antibody binding for the immobilized receptor, reported as the number of antibody molecules (Ab) bound per EcAII molecule (antigen; Ag) with Kidrolase (black), N21-EcAII (green), N26-EcAII (blue), and C8-EcAII (red). Data are presented as the mean values, and the error bars indicate standard deviation.

Sensing efficiency

The sensing efficiency is defined as the number of analyte molecules detected per molecule of the immobilized receptors or Ab/Ag binding ratio. It was observed that the native tetrameric EcAll has a significantly better sensing efficiency than the Histagged variants despite having low surface density and poor sensitivity. The average of 2.3 target antibody molecules were detected per immobilized Kidrolase molecule with binding ratio for about 2 for 1. This is 2 to 3-fold more efficient than for the tetrameric Histagged EcAll variants, which can only bind to 0.9 antibody per molecule. This is explained by the higher density on the sensing surface than Kidrolase (Table 1). The average density of immobilized Kidrolase was ~7fold (1.6 x 10¹¹ molecules cm⁻²) lower than native EcAII with average surface density of ~ 1.15×10^{12} molecules cm⁻² in the crystal lattice (PDB 3ECA) plane. This gives an average distance between protein tetramers of 11.4 nm, from center to center (c.t.c) in comparison to average c.t.c distance between the immobilized tetramers of 2.7-fold greater (30.5 nm). On the contrary, the surface density tetrameric Hist-tagged EcAll variants immobilized on the sensing surface was

1.3 to 1.8 fold greater than in crystal lattice, with c.t.c distances between tetramers of about ~1.2 to 1.3-fold

shorter. This could be caused by the higher packing with the 8- and 25-residue terminal linked tetramers to favor bivalent antibody binding. Nonetheless, the lower sensing efficiency of the His-tagged EcAII is compensated by the greater surface density and by the reproducibility which ultimately ensure greater sensitivity when detecting anti-EcAII in serum.

The P4SPR advantage

Affinité instruments' P4SPR is a modular, portable, multichannel system with flexibility in assay development and optimization. Its accessibility and ease of use allow for quick evaluation of assay parameters such as the effect of immobilization approaches on receptor surface coverage as well as sensing efficiency. Its multichannel feature increases result precision and additionally, the demonstrated assay sensitivity proves the system to be a powerful analytical tool for immunosensing which is especially interesting for analysis of biotherapeutic agents in biological samples. No required sample pre-treatment and minimal and shorter testing procedure are one of the main advantages over conventional ELISA technique.

Advantages of P4SPR compared to ELISA

- Label-free
- Minimal reagents used
- Fewer preparation steps
- Reduced time and reagent costs

About Affinité Instruments

Established in 2015 as a spin-off of the Université de Montréal, Affinité instruments' foundation is built deep knowledge accrued on over a decade of research in SPR. The commercialization of promising innovations is spearheaded by a diverse leadership experienced in business, science and engineering.

Additional information

This is an open access article published under an ACS AuthorChoice License, which permits copying and



redistribution of the article or any adaptations for noncommercial purposes.

Annex - Surface coverage of immobilized EcAII

The surface coverage (Γ , ng cm⁻²) of the immobilized native or His-tagged EcAII was calculated from the change in the wavelength ($\Delta\lambda$ SPR) upon immobilization of EcAII using the following equation:

$$\Gamma = \rho(-l_d/2) \ln(1 - (\Delta \lambda / m(n_{\text{SAM}} - n_{\text{medium}})))$$

where ρ is to the density of the adsorbed protein monolayer (1.3 g cm⁻³), I_d is the plasmon penetration distance (~230 nm), $\Delta\lambda$ is the shift in the wavelength associated with protein immobilization, m is the refractive index sensitivity of the SPR sensor (1765 nm/RIU), n_{SAM} is the refractive index of the peptide SAM (1.57 RIU), and nmedium is the refractive index of the buffer (1.33476 RIU). The total amount of immobilized EcAII on the sensing surface (Q) was determined using the formula Q = Γ S, where S = 0.166 cm² in contact with the protein.

where p is to the density of the adsorbed protein monolayer (1.3 g cm–3), ld is the plasmon penetration distance (~230 nm), $\Delta\lambda$ is the shift in the wavelength associated with protein immobilization, m is the refractive index sensitivity of the SPR sensor (1765 nm/RIU), n_{SAM} is the refractive index of the peptide SAM (1.57 RIU), and n_{medium} is the refractive index of the buffer (1.33476 RIU). The total amount of immobilized EcAII on the sensing surface (Q) was determined using the formula Q = Γ S, where S = 0.166 cm² in contact with the protein.